













Hamilton T1 Ventilator Standard Operating Procedure

Adult (>15Kgs) Invasive Ventilation Set-up and Operation

Adapted from SJH ED SOP: Author Ross Archibald (NHSL employee)

Instruction		Example / Note
1.	<p>Items required:</p> <p>Circuit pack (code: 260207) containing;</p> <ul style="list-style-type: none"> - Disposable coaxial circuit - Disposable proximal flow sensor - Disposable expiratory valve 	
2.	<p>a) Connect AC power cable to wall and turn on</p> <p>b) Connect oxygen hose to cylinder or wall outlet (if supplementary oxygen is required).</p> <p>c) Turn power ON (hold power button briefly and allow auto start up procedure to complete).</p> <p>d) The ventilator will run a short self-check then displays the standby screen. Ensure no major faults indicated (RED flashing light and accompanying error message).</p>	
3.	<p>a) Open circuit pack.</p> <p>b) Detach expiratory valve from disposable circuit. <u>Ensure rubber membrane cap remains attached to expiratory valve.</u></p> <p>c) Insert expiratory valve into expiratory port (membrane cap first) and rotate clockwise until click felt.</p>	
4.	<p>Connect coaxial circuit ends to inspiratory port and expiratory port as shown.</p>	
5.	<p>Separate ends of proximal flow sensor tubing as shown.</p>	
6.	<p>Connect proximal flow sensor tubing to ventilator (blue to blue, clear to white).</p>	

7.	<p>Select 'Preop check' on standby screen (bottom left). Perform the leak and flow sensor checks.</p>	
8.	<p><u>Leak test</u></p> <p>a) Don gloves.</p> <p>b) Remove red cap from end of circuit.</p> <p>c) Select "Leak test" and follow instructions on screen.</p> <p>d) Block end of flow sensor (dark blue part) and release when instructed.</p> <p>GREEN tick indicates successful test. RED cross requires repetition of above steps.</p>	
9.	<p><u>Flow sensor check</u></p> <p>a) Select "Flow Sensor" and follow instructions on screen.</p> <p>b) When prompted, attach calibration adapter (opaque part) to the end of the flow sensor.</p> <p>c) When prompted, flip both parts 180° so the calibration adapter is directly connected to the coaxial circuit.</p> <p>d) Remove calibration adapter (don't throw away) and return flow sensor adaptor to original position as shown.</p> <p>GREEN tick indicates successful test. RED cross requires repetition of above steps.</p>	
10.	<p><u>O2 sensor check</u></p> <p>Not routinely required unless indicated by ventilator. Perform check if not done within the last 4 weeks.</p> <p><u>Loss of power check</u></p> <p>Perform a loss of power check by disconnecting the ventilator from the mains power supply. Verify that the alarm is activated and that the ventilator is powered by its backup battery. Reconnect the ventilator to the AC power supply</p>	
12.	<p>Once three GREEN ticks appear on Pre-op check menu*:</p> <p>Connect ventilator tubing as per departmental protocol (ETT > catheter mount > HME Filter > capnography > ventilator circuit).</p> <p>Select desired ventilator settings and input patient specific information (sex, height).</p>	
<p>For further information, support or queries please contact: Jake Maxwell, UK Manager, Eden Medical UK 07551 171280 jake.maxwell@edenmedical.co.uk</p>		

***If the pre-operational checks fail**

1. Check breathing circuit – no disconnection, properly connected and no leakage
2. Check flow sensor – properly connected (as per steps 5-6)
3. Check expiratory valve is properly installed and seated (step 3)
4. Replace breathing circuit, flow sensor and expiratory valve set – ensure properly connected and seated
5. Return ventilator to medical physics for servicing