Index Number DCN/MRI/SOP/	Procedure for scanning (ITU) patients in DCN MF		al care	NHS Lothian
Site	Dept of Clinical Neurosciences	Modality	MRI	
Department	DCN X-Ray	Location	DCN MRI	

# **Department of Clinical Neurosciences** Neurological MRI under General Anaesthesia **SOP for Critical Care**

Produced by L Steel, R Sutherland, S Jamieson, M Schwarz

Significant injury and death can occur in the MRI as a consequence of poor understanding of the risks of MRI. The new 3 Tesla MRI in DCN has a magnetic field twice as powerful as our previous MRI, providing extremely high quality imaging.

This SOP is designed for you and your patients' safety, please complete every step carefully. ALL STAFF REQUIRED TO WORK IN MRI MUST HAVE RECEIVED MRI SAFETY TRAINING

## Safety Assessment for MRI

'The DCN GA MRI test	<b>'</b> – Lie your patient <i>compl</i>	letely flat i	for 90 minutes
----------------------	--	---------------	----------------

'The DCN GA	A MRI test' – Lie your patient <i>completely flat</i> for 90 minutes		
	Has there been a significant desaturation?	Yes	No 🔙
•	Has the FiO2 increased or higher than 0.5?	Yes	No
•	Has the cardiovascular status become unstable at any point?	Yes	No 🗌
•	Does your patient have a chest drain that cannot be clamped?	Yes	No 🗌
•	Is there a high vasopressor/ionotrope requirement?	Yes	No 🗌
If the	answer to any of the above is yes, your patient is not currently fit	for MRI	
<ul><li>YOUF</li></ul>	R PATIENT MAY NOT BE IMMEDIATELY ACCESSIBLE FOR UP TO 90	MINUTES	
How to B	ook a Neurology MRI:		
■ R	equest the scan on TRAK – MRI Brain or Spine only please		
	<ul> <li>Please avoid using N/A in the safety questions – this will res being unable to authorise the scan</li> </ul>	ult in the Radiolo	ogist
• D	iscuss the scan with the on call Consultant Neuroradiologist (501	.19 Or via Switch	)
• C	nce authorised, discuss MRI slot availability with the <b>DCN MRI Ra</b>	diographers (508	370)
• C	all the <b>on-call Neuroanaesthetist</b> to organise Anaesthetic cover (8	519 via DCN pag	er)
	n mind there may not be capacity for same-day scanning, but we v e the scan as soon as possible	vill endeavour to	

Date Issued: Version No. Author: 20/09/2020 1.0 Stuart Jamieson	Authorised by: Stuart Jamieson	Date of Last Review: 11/02/2022	Next Review Date: 11/02/2024	Page 1
---	--------------------------------------	---------------------------------------	------------------------------------	--------

Index Number	Procedure for scanning ventilated critical care	NHS Lothian
DCN/MRI/SOP/	(ITU) patients in DCN MRI	

### Team for MRI Transfer

- MRI-trained Anaesthetist
- MRI-trained ODP
- ICU nurse
- Radiography team
- All staff must complete their own MRI safety questionnaires before entering MRI
- Please avoid extraneous staff as this is a significant safety issue

## Preparation for Transfer to MRI Checklist

•	MRI safety questionnaire complete for the patient	
•	Roll or log roll to check for metal under the patient	
•	Remove all metal items (hair clips, jewellery etc)	
•	Reduce all equipment to the minimum necessary	
•	Reduce number of IV infusions to minimum necessary	
	Replace all IV lines with extra long extensions	

**If you are unsure about any item**, please discuss with MRI radiographers (ex 50870 / 50906) *before* transfer to MRI

## Arrival at DCN MRI (1st Floor) (Intraoperative scanner in DCN theatre)

- Ring the bell to alert the MRI Radiographers that you have arrived
- In the bay outside MRI transfer to the MRI table
- Remove all MRI Unsafe equipment and monitoring
- Attach MRI safe monitoring
- Detach portable ventilator and connect c-circuit
- Enter the MRI prep room patient head first next to the scanner can reconnect to MRI ventilator
  - o Remove all metal/ferrous items from your person and leave it in this room
  - Disconnect ICU pumps from patient
  - o Complete the GA MRI Safety Pause just prior to transferring into scan room
- Separate drug pumps should be set up in the control room, ready to be connected to the extra long IV lines.
- Leave the original ICU pumps on the patient ICU bed outside the MRI suite
  - o All IV pumps are MRI Unsafe and CANNOT ENTER THE SCANNING ROOM

Date Issued: Version No. 20/09/2020 1.0	Author: Authoris Stuart Stuart Jamieson Jamiesor	Review:	Next Review Date: 11/02/2024	Page 2
---	--	---------	------------------------------------	--------

### Transferring the Patient into the Scan Room

- Once MRI safety pause complete can transfer into scan room patient feet first
- Connect to MRI ventilator in scan room
- Connect IV lines to lines already fed through waveguide
- Radiographers will raise scanner table to lift patient off MRI trolley up to scanning position
  - o Ensure no lines / tubing etc get caught
- Hearing protection given to patient
- Arm supports and straps should be used to prevent lines getting caught
- For brain imaging ventilator tubing may need adjusting to enable MRI head coil to be positioned

A **test run** of moving patient through the magnet should be carried out to check lines / tubing will not get pulled.

#### After the scan

- Once scan is finished patient can be moved direct from the scanning room to the ICU bed outside MRI
- Reconnect ICU pumps, monitoring and portable ventilator outside MRI
- ODP / Radiographers remember to recharge MRI monitoring batteries and plug in mobile monitoring unit and check if Oxygen cylinders needing replaced.

### Scanning on the ground floor DCN MRI scanners

Occasionally ventilated scans will need to happen on the ground floor scanners if the intraoperative scanner is not available or a 1.5T scan is needed. Preference is always for the intraoperative scanner.

- As the ground floor environment will be less familiar to anaesthetic and ODP staff, when a
  case is arranged to happen on the ground floor the anaesthetist and ODP should visit the
  area to understand the environment and patient pathway.
- Patient arrives to DCN Imaging inpatient holding bay and is transferred onto the MRI table
   (Radiographer to ensure correct table for scanner to be used)

Date Issued: Version No. 20/09/2020 1.0	Author: Stuart Jamieson	Authorised by: Stuart Jamieson	Date of Last Review: 11/02/2022	Next Review Date: 11/02/2024	Page 3
---	-------------------------------	--------------------------------------	---------------------------------------	------------------------------------	--------

- Remove all MRI Unsafe equipment and monitoring
- Attach MRI safe monitoring
- Detach portable ventilator and connect c-circuit
- Move patient from the inpatient holding bay round into the MRI controlled access are to the MRI prep room – patient head first – can reconnect to MRI ventilator
  - o Remove all metal/ferrous items from your person and leave it in this room
  - o Disconnect ICU pumps from patient
  - o Complete the GA MRI Safety Pause just prior to transferring into scan room
- Separate drug pumps should be set up in the control room, ready to be connected to the extra-long IV lines.
- Move the original ICU pumps back to the patient ICU bed in the inpatient holding area
  - o All IV pumps are MRI Unsafe and CANNOT ENTER THE SCANNING ROOM
- Patient taken into scan room on successful completion of GA MRI Safety pause
  - Transfer process on ground floor is same as described for first floor intraoperative scanner.

#### Other Points to Consider

- For any emergency the patient must be removed from the scan room and ideally outside the MRI controlled access area.
- The MRI scanner is always on and so there is always a strong magnetic field nearby
- Some types of tattoos or permanent makeup contain ferromagnetic material
- Please bear in mind the maximum weight capacity 250Kg for the MRI when booking a scan. The diameter of the bore is 70cm if there is doubt whether the patient will fit in MRI please contact the Radiographers (50870) who can come and assess.
- Moving within a powerful magnetic field can cause staff to experience a metallic taste, headaches or sensory disturbance when very close to opening of magnetic bore – please alert someone if you feel unwell.
- MRI infusion pumps these are now available for use to those trained. Remember to remove the black battery from the mobile unit on the drip stand and replace in the charger after use. Battery should not be removed in the scanning room.

Date Issued: Version No.   Author:   Authorised by:   Date of Last   Next Review   Page 20/09/2020   1.0   Stuart   Stuart   Stuart   Stuart   Tamieson   Stuart   Stuart   Tamieson   11/02/2022   11/02/2024   Stuart   Tamieson   Stuart   Tamieson   Stuart   Tamieson   Stuart   Tamieson   Stuart   Tamieson   Stuart   Stuart	Date Issued: 20/09/2020
---	----------------------------