

LUHT Critical Care Directorate Guidelines

TED Anti-Embolism Stockings SIZING GUIDE:

Below Knee:

1. Measure around widest part of calf

Small	< 33 cm
Medium	33-38 cm
Large	38-40cm
Extra Large	> 40cm

Thigh Length:

1. Measure calf circumference to determine the appropriate size
2. Measure the upper thigh circumference BUT DO NOT APPLY thigh length stockings if circumference exceeds the following:

Small	> 65cm
Medium	> 70cm
Large	> 78cm
Extra Large	> 82cm